

# Keysight Technologies 8500 FE-SEM

Compact System for Low-Voltage,  
High-Performance Imaging



Data Sheet

## Overview

The new Keysight Technologies, Inc. 8500 FE-SEM offers researchers a field emission scanning electron microscope (FE-SEM) right in their own laboratory. This compact, innovative system has been optimized for low-voltage imaging, extremely high surface contrast, and resolution typically found only in much larger and more expensive field emission microscopes.

The 8500 is easy to install and use. No dedicated facilities are required, only an AC power outlet. About the size of a laser printer, the scientific-grade system provides capabilities to researchers in ordinary labs that previously were only available with conventional FE-SEMs installed in centralized facilities. The 8500 has been engineered to deliver consistent, reproducible performance and the industry's lowest total cost of ownership for a FE-SEM.

## Enhanced Imaging Capabilities

The 8500 FE-SEM offers several imaging techniques for enhancing surface contrast and allowing nanoscale features to be observed on a wide variety of nanostructured materials, including polymers, thin films, biomaterials, and other energy-sensitive samples on any substrate, even glass.

The system's continuously variable imaging voltage is tunable from 500 to 2000 volts as an operational parameter, rather than a setup choice. The 8500 eliminates charging of nonconductive samples without the need to coat the samples, which can mask nanoscale features, or resort to increased pressure operation, which can degrade resolution.

In addition, the 8500 utilizes a four-segment microchannel plate (MCP) detector that provides topographic imaging along two orthogonal directions to enhance surface detail. This technique has been demonstrated to clearly resolve sub-nanometer atomic steps on the surface of crystalline substances such as polytype 6H-SiC

## Novel Design

The core technology inside a scanning electron microscope is the electron beam column, which extracts, collimates, shapes, scans, and focuses the electron beam. A conventional electron beam column relies on a combination of precision-machined electromagnetic and electrostatic elements to control the electron beam. The coils that form the critical elements are typically hand-wrapped by technicians to achieve uniform electromagnetic fields. Moreover, closed-loop cooling and sophisticated vibration isolation are often required at the system level to manage the high currents in the lenses and other elements. The resultant high-resolution electron beam column is both large and expensive.


Silicon-based microfabrication techniques enable Keysight Technologies to design and fabricate a miniature electrostatic electron beam column combined with a thermal field emission electron source. The 8500 FE-SEM design utilizes stacks of silicon on insulator to form all of the lenses, apertures, and deflectors in the electron beam column. Patented technology allows these electron beam columns to be built wafer-scale on

## Features and Benefits

- Resolution and imaging equal to that of conventional FE-SEMs
- Variable low voltage eliminates charging and the need for sample coating
- Programmable X, Y, Z stage allows user to set precise coordinates, scan, and save information
- Miniature electrostatic lens design ensures repeatable performance without constant re-tuning
- Compact size enables easy installation in any research laboratory and does not require special facilities

## Applications

- Polymers
- Thin films
- Biomaterials
- Nonconductive samples
- Energy-sensitive materials
- Glass substrates



Schematic of electrostatic electron beam column.

150mm substrates. The columns are fabricated with the precise aperture diameters and repeatable alignment tolerances required to minimize aberrations that can degrade image quality.

The system's thermal field emission electron source, meanwhile, provides high brightness, high stability, small virtual source size, and low energy spread and long-lasting consistent performance. A quad-segmented MCP detector is located just below the objective lens of the electron beam column directly and above the sample. This detector collects both backscattered and secondary electrons. The MCP may be operated either in a standard mode, where all the channels are added together, or in a differential mode (topographic mode), in which opposite sides of the detector are dynamically subtracted.

## Performance Advantages

As described above, the 8500 FE-SEM embodies the successful miniaturization of the core technology in a scanning electron microscope. The 8500 is optimized for low-voltage imaging and sub-10nm resolution. Its thermal field emission electron source provides high signal-to-noise ratios and consistent, long-lasting performance, while secondary and backscatter electron detection capabilities provide a rich data set for each sample.

Furthermore, the electrostatic lens design of the 8500 delivers repeatable performance without the constant re-tuning necessitated by the hysteresis in magnetic lenses found in conventional SEMs. This design allows researchers to store and return to any operating setup immediately with negligible fine adjustment, making the 8500 an ideal choice for multiple-user environments.

## Ultimate Ease of Use

Just as the factory-calibrated 8500 FE-SEM is approximately the same size as a laser printer, it also offers similar plug-and-play performance. In fact, the compact size of the 8500 facilitates easy installation in practically any existing lab that has an AC power outlet.

Sample preparation and loading is both simple and fast. An X, Y, Z programmable stage lets users set specific coordinates, scan and then store the locations to repeat experiments with precision and confidence. For optimum control, the system's powerful software package features an intuitive graphical user interface (GUI) designed for novice and expert users.

The highly versatile 8500 can speed research and product development cycles by decreasing the time from experimentation to characterization. With the 8500, convenient access to FE-SEM performance is now a reality.

## ECD Cartridge

The electron source, the electron beam column, and the electron detector of the 8500 FE-SEM are all combined in a field-replaceable ECD cartridge. When the long-lasting electron source is finally depleted, the entire ECD cartridge can be replaced on-site to provide the 8500 with not only a new source of electrons, but a new pre-aligned electron beam column and a new MCP detector – essentially a brand new FE-SEM.

Unlike conventional FE-SEMs, the electron source of the 8500 is easy to turn off and start up again quickly for imaging. By powering off the electron source when it is not in use, the lifetime of the source can be further extended, thus lowering the system operating costs and delivering the industry's lowest total cost of ownership for a FE-SEM.




Figure 1. Bacillus megaterium incubated with L-omwaprin.




Figure 2. High resolution image of gold islands on carbon.



Figure 3. Low voltage image of uncoated electrode with CNTs on tip.




Figure 4. High resolution image of celgard.

## 8500 FE-SEM Specifications

### Performance

Beam voltage:	500 to 2000V
Beam current:	0.2 to 1nA
Resolution:	<10nm at 1000V
Magnification*:	250 to 65,000 X
Digital zoom:	10X
Scan field:	1 x 1mm (max)
Electron source:	Schottky field emission
Detector modes:	SE, BSE, Topo

### Sample

Sample size:	100 x 60mm (max)
Sample thickness:	30mm (max)
Viewable area:	50 x 30mm (max)
Sample mounts:	Standard SEM stubs
Electrical activation:	Vacuum feedthrough

### System Control

PC:	Windows 7
Basic user interface:	Simple image capture and controls
Expert user interface:	Expert image optimization

### Image

Image formats:	JPEG, TIFF, BMP, PNG
Image resolution:	User selectable up to 2048 x 2048 pixels
Scan rate:	Slow scan to video rate
Noise reduction:	Frame and pixel averaging

### Motorized Stage

Piezo stage:	1 µm accuracy
X, Y, Z travel:	50 x 50 x 10mm

### Vacuum System

Chamber vacuum:	1 e-4 Torr
Pumpdown time:	3 minutes
Turbo pump:	80 liters per second
UHV pump:	Ion pump with gettering

### Dimensions

Microscope:	584(W) x 470(D) x 584(H) mm; 72 kg
Pump unit:	203(W) x 254(D) x 203(H) mm; 4 kg

### Installation Requirements

Power:	100/120/220–240 VAC; 50/60 Hz
Operating temperature:	5 to 40°C
Humidity:	20 to 80% RH
Compressed air:	Not required
Dry nitrogen:	Not required**
Water cooling:	Not required

\* Electron optical magnification relative to a 3nm pixel on a VGA monitor.

\*\* A N2 vent port is configured and available.

## Nano Mechanical Systems from Keysight Technologies

Keysight Technologies offers high precision, modular nano-measurement solutions for research, industry, and education. Exceptional worldwide support is provided by experienced application scientists and technical service personnel. Keysight's leading-edge R&D laboratories ensure the continued timely introduction and optimization of innovative, easy-to-use nanomechanical system technologies.

[www.keysight.com/find/nanoindenter](http://www.keysight.com/find/nanoindenter)

For more information on Keysight Technologies' products, applications or services, please contact your local Keysight office. The complete list is available at:

[www.keysight.com/find/contactus](http://www.keysight.com/find/contactus)

### Americas

Canada	(877) 894 4414
Brazil	55 11 3351 7010
Mexico	001 800 254 2440
United States	(800) 829 4444

### Asia Pacific

Australia	1 800 629 485
China	800 810 0189
Hong Kong	800 938 693
India	1 800 112 929
Japan	0120 (421) 345
Korea	080 769 0800
Malaysia	1 800 888 848
Singapore	1 800 375 8100
Taiwan	0800 047 866
Other AP Countries	(65) 6375 8100

### Europe & Middle East

Austria	0800 001122
Belgium	0800 58580
Finland	0800 523252
France	0805 980333
Germany	0800 6270999
Ireland	1800 832700
Israel	1 809 343051
Italy	800 599100
Luxembourg	+32 800 58580
Netherlands	0800 0233200
Russia	8800 5009286
Spain	0800 000154
Sweden	0200 882255
Switzerland	0800 805353
	Opt. 1 (DE)
	Opt. 2 (FR)
	Opt. 3 (IT)
United Kingdom	0800 0260637

For other unlisted countries:

[www.keysight.com/find/contactus](http://www.keysight.com/find/contactus)

This information is subject to change without notice.

© Keysight Technologies, 2014  
Published in USA, November 13, 2014  
5990-6183EN

[www.keysight.com](http://www.keysight.com)